









# **ENDODERMIC MASSAGE** (VACUUM+ROLLERS) + R.F. + SOFTLASER



- 2 BODY HANDPIECES with built-in display and controls for easy adjustment and display of treatment parameters. Excellent ergonomics to facilitate the job of the beautician.
- · Improved skin elasticity and appearance.
- · Body Slimming.
- Firming.
- · Targets blemishes due to fat and cellulite.
- 1 FACE HANDPIECE with medium applicator for anti-ageing treatments and skin rejuvenation.
- 1 FACE HANDPIECE with small applicator dedicated to reducing wrinkles and treating stretch marks.









**RADIOFREQUENCY** non-invasive treatment that uses the action of electromagnetic waves to counter various types of imperfections including: skin laxity, stretch marks and even improvement of imperfections related to cellulite.

#### THE ENDODERMAL TREATMENT

is a vacuum connective tissue massage, i.e. it aims to reactivate the circulation and oxygenate the cells by means of a mechanical massage that stresses the skin, creating a vacuum and massaging it with metal rollers.

THE SOFT LASER acts on the skin by emitting light at a wavelength of 940nm. The laser energy reaches the subcutaneous adipose tissue and stimulates it.













## software with the possibility of customising treatments.

• SAFETY AND SANITATION

Each handpiece features removable and sanitised applicators.

INCLUDING NO.2 DRY OIL WITH INNOVATIVE AND SPECIFIC BODY AND FACE FORMULAS. SUITABLE FOR TREATMENT WITH ENDODERMIC MASSAGE EQUIPMENT.

IT ALLOWS THE HANDPIECE TO GLIDE MORE SMOOTHLY, IT NOURISHES AND MOISTURISES THE SKIN AND DOES NOT LEAVE IT GREASY.

N. 1 ITEM 69C - Dry Oil with Phosphatidylcholine Body - 500ml.

Specific for treatment of skin imperfections in people with localised fat deposits and/or cellulite.

 $\ensuremath{\mathsf{N.\,1}}$  ITEM 208V - Dry Oil with Multivitamin Synergy for Treatment

Face - 200ml

Specific for anti-ageing treatments to prevent skin ageing, protects and tones the skin, preserving it from laxity and wrinkles.







ART. 9157

#### **TECHNICAL FEATURES**

- Vacuum from 10 to 80kPa
- from 0 to 60 joule/cm<sup>2</sup>
- Soft Laser from 0 to 20W
- Soft laser wave length 940nm
- Roller Speed from 0 to 36 rpm

#### **TECHNICAL FEATURES**

- Power supply 220 - 240V ~ 50/60Hz - 750W
- Device size 64x61,5 H108,5cm
- · Device size with handpiece holder mounted 64x61,5 H168cm
- · Package sizes 75x62x123cm
- · Device weight complete with handpieces 54kg
- Device weight with packaging 84kg



# **EXPO MÜSTER AREA SUD**

Via Massafra, 30i 74015 Martina Franca (TA) - Italy Tel. 080-4838400 e-mail: expo.sud@muster-dikson.com

# **♦**MÜSTER & DIKSON S.p.A.

## SEDE, STABILIMENTI ED ESPOSIZIONE

Via Privata da via Kennedy snc 20023 Cerro Maggiore (MI) - Italy Tel. 0331-525111 (r.a.) e-mail: sales@muster-dikson.com

Export Department: Tel. +39 0331-525111 (r.a.) e-mail: export@muster-dikson.com

www.muster-dikson.com

# MÜSTER & DIKSON ESPAÑA

Avd Gran Vía de les Corts Catalanes n° 524, ent 2° 08011 BARCELONA - España Tel. +34 93 421 0589 e-mail: carmen.gonzalez@muster-dikson.com

#### MÜSTER & DIKSON NORTH AMERICA INC.

350 Harry Walker Parkway North - Suite 14 Newmarket Ontario - L3Y-8L3 Phone +1 905.235.5580 - Fax +1 905.235.5524 e-mail: info@dikson.ca - www.dikson.ca









